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电子邮件：service@genepharma.com

©2005 - 2006 Genepharma Corporation. All rights reserved. For research use only. Not intended for any animal or human therapeutic or diagnostic use.

Kit Contents and Storage

Cat. No. QPM-050	Size: 50rxns
Cat. No. QPM-051	Size: 100 rxns

U6 snRNA Real-time PCR Normalization Kit Reagents

U6 snRNA Real-time PCR Normalization Kit Reagents box includes the following items. Store the components at -20°C . **Reagents must be stored in the dark.**

Reagent	Amount		Storage
	50rxns	100rxns	
Real-time PCR Master Mix ($2\times$) *	1.0ml	2.0 ml	4°C
U6 snRNA primer Set ($5\mu\text{M}$)	$80\mu\text{l}$	$160\mu\text{l}$	-20°C
U6 snRNA RT primer ($10\mu\text{M}$)	$15\mu\text{l}$	$30\mu\text{l}$	-20°C
Taq DNA polymerase ($5\text{U}/\mu\text{l}$)	$20\mu\text{l}$	$40\mu\text{l}$	-20°C
Sterilized ddH ₂ O	1.0 ml	1.0 ml	4°C

*Include dNTP Mixture, Mg^{2+} , SYBR Green I, Rox.

Accessory Products

Accessory Products

Some of the reagents supplied in the U6 snRNA Real-time PCR Normalization Kit as well as other products suitable for use with the kit are available separately from Genepharm. Ordering information is provided below. For more information, refer to our Web site (www.genepharm.com).

Item	Amount	Catalog no.
Hairpin-it TM miRNAs qPCR Quantitation Kit	$20\mu\text{l} \times 50\text{ rxns}$	QPM-010
	$20\mu\text{l} \times 100\text{ rxns}$	QPM-011

Introduction

Overview

Introduction

The U6 snRNA Real-time PCR Normalization Kit is a sensitive and specific method using real-time PCR for the normalization and profiling of microRNAs (miRNA) from total RNA samples.

MiRNAs are small, single-stranded, ~19–23 nt RNA molecules encoded in the genomes of plants, animals, and viruses. Mature miRNAs enter the RNA-induced silencing complex (RISC) and guide the RISC to induce translational repression or endonucleolytic cleavage of specific target mRNAs.

The U6 small nuclear RNA (snRNA) undergoes major conformational changes during the assembly of the spliceosome and catalysis of mRNA splicing. It associates with the specific protein Prp24p, and a set of seven LSm2p-8p proteins, to form the U6 small nuclear ribonucleoprotein (snRNP). The spliceosome performs pre-mRNA splicing in eukaryotes.

Description of U6 snRNA Real-time PCR Normalization Kit

The Real-time PCR normalization assay of miRNAs using U6 snRNA as a housekeeping gene includes two steps, U6 snRNA RT reaction and real-time PCR. The U6 specific RT primer is supplied in the U6 snRNA Real-time PCR Normalization Kit, the RT-PCR kit (Cat # QPG-090) is available separately.

Users can use the U6 snRNA Real-time PCR Normalization Kit to normalize the expression ratio of a specific microRNA relative to the U6 snRNA between different samples.

Basic of Relative Quantification Calculation

Relative quantification of microRNA gene expression is calculated from the threshold cycle (C_T) values obtained in experiments. The threshold cycle is the cycle at which a statistically significant increase in ΔRn is first detected. Thus, wells with higher initial template concentrations reach the threshold value at lower cycle numbers during PCR than wells containing lower initial template concentrations.

For an ideal microRNA quantification assay, each cycle in the PCR reaction corresponds to a twofold increase in PCR product. Therefore a change in threshold cycle number of 1 equates to a two fold difference in initial template concentration. For example, if the C_T increases by one, the initial template concentration decreases two fold. Conversely, if the C_T decreases by one, the initial template concentration doubles. This property is used in the calculation of relative quantification values for the cytokine and internal control target sequences in each well.

For more information, see the **Data Analyzing** chapter in this manual.

Components of the Kit

Introduction

This section provides more information about the reagents supplied in the U6 snRNA Real-time PCR Normalization Kit.

Real-time PCR Master Mix 2× Conc.

The 2× Real-time PCR Master Mix consists of a proprietary buffer system, MgCl₂, dNTPs, SYBR Green I and Rox reference dye. The mix has been confirmed to be work well in restrict lab experiments.

IMPORTANT

It's very important for you to store the 2× Real-time PCR Master Mix in the dark.

U6 snRNA primer set 5μM

The pre-designed U6 snRNA primer set has been confirmed to be work well in restrict lab experiments. Data base information for U6 snRNA is provided below.

<u>Species</u>	<u>GeneBank Accession#</u>	<u>Primer Tm</u>	<u>Amplicon length</u>
<i>Homo sapiens</i>	M14486	58~62℃	70~100 bp

U6 snRNA RT primer 10μM

The pre-designed U6 snRNA RT primer is provided as a 10 μM solution. We recommended the final concentration of the RT primer 60nM for a reverse transcription reaction. **Please dilute the RT primer for 10 folds first and use 1.2μl per RT reaction.**

Taq DNA polymerase 5U/μl

For most PCR reaction, the final concentration of DNA polymerase was usually 0.05U/ μ l per reaction volume.

User-Supplied Required Materials

Reverse Transcriptase

For detection and quantification of the miRNAs in your sample, a RT reaction use the miRNA specific RT primer provided in the kit was needed. We recommended using traditional Reverse Transcriptase such as MMLV.

RNase Inhibitor

To keep the activity of Reverse Transcriptase, adding RNase Inhibitor in your RT reaction was recommended. The final concentration usually was 0.25 U/ μ l.

Rox reference dye

To normalize the fluorescence difference from well to well, adding Rox reference dye in your qPCR system was recommended when using the Real-Time PCR instruments such as ABI PRISM 7000/7300/7500/7900, MX3000p/4000.

Optical PCR plates or tubes

As to eliminate the impact on fluorescence detection caused by PCR tubes or plates, using optical PCR plates and tubes was recommended to corporate the requirements of the Real-Time PCR instruments such as ABI PRISM 7000/7300/7500/7900, MX3000p/4000.

Method

Performing U6 snRNA RT Reaction

Introduction

The Hairpin-itTM miRNAs real-time Quantitation assay of miRNAs includes two steps, stem-loop RT reaction and real-time PCR detection. To profiling microRNA's expression relative to a said housekeeping gene, U6 snRNA or 5S rRNA qPCR Normalization Kit should be needed. The U6 snRNA normalization kit is used to normalize the expression level of microRNAs relative to U6 snRNA, which is a small nuclei RNA abounded in cells. Same volume RNA sample of microRNA assay and U6 normalization assay is asked to add in the RT reaction system, then, the RT product is quantified using real-time PCR.

Handling U6 snRNA RT Primer

Please dilute the RT primer for 10 folds first and use 1.2μl per RT reaction.

Preparing The U6 snRNA RT Reaction Mix

The template for the U6 snRNA Real-time PCR Normalization Kit can be total RNA or cell lysate. The RNA template's input should be at least 1μg or more as for as the requirement of the experiment. Other than the stability at 42°C, there is no special requirement of Reverse Transcriptase for the RT reaction. We recommended a standard 20μl RT reaction size.

The following table provides RT reaction Mix volumes for a standard 25 μl reaction size.

Component	Final Con.	Vol /1 rxns
5×RT Buffer	1×	4μl
dNTP (10mM)	0.375mM	0.75μl
U6 snRNA/miRNA RT Primer (1μM)	60nM	1.2μl
RNasin (40U/μl)	0.5U/μl	0.25μl
MMLVReverse Transcriptase (200U/μl)	40U	0.2μl
RNA Sample (Total RNA) *	1μg	Xμl
RNase Free H ₂ O	—	To 20μl

* It's important to keep the RNA sample volume or quantity of U6 normalization RT assay and miRNA RT assay coincidence. Mix the RT reaction reagents and Mix before RT reaction.

Continued on next page

Performing U6 snRNA RT Reaction, continued

Performing miRNAs RT Reaction

Standard RT Reaction Program

30 min at 16_°C, 30~45 min at 42_°C, 5 min at 85_°C

NOTE

Keep all components, reaction mixes and samples on ice. After assembly, transfer the reaction mixes to a thermal cycler preheated to the cDNA synthesis temperature and immediately begins RT reaction.

Performing miRNAs & U6 snRNA Real-Time PCR Reaction

Handling the U6 snRNA RT product

For the U6 snRNA Real-time PCR Normalization Kit, pipet 2 µl U6 snRNA RT reaction product as the template for real-time PCR step subsequently. Store the surplus microRNA RT reaction product at -20°C.

Preparing Real-Time Reaction Mix

The following table provides Real-time PCR reaction Mix volumes for a 20 µl reaction size. Note that three repeat tubes be performed in an experiment at least.

Component	Final Con.	Vol /1 rxns
2×Real-time PCR Master Mix ¹	1×	10 µl
U6 snRNA Primer set(5 µ M) ²	0.2 µ M	0.8 µl
U6 snRNA RT product		2 µl
Taq DNA polymerase (5U/µl)	0.5 U/µl	0.2 µl
dd H ₂ O		To 20 µl

¹The 2×Master Mix contains Mg²⁺, dNTP,SYBR Green I and Rox dye.

²The U6 snRNA Primer set includes PCR forward and reverse primers.

For multiple reactions, prepare a master mix of common components, add the appropriate volume to each tube or plate well, and then add the unique reaction components.

Performing miRNAs & U6 snRNA Real-time PCR Reaction

Program the real-time PCR instrument to perform PCR amplification as shown below.

95°C for 3minute hold, 40 cycles of:

95°C, 15 seconds ,62°C, 40~60seconds.

Continued on next page

**NOTE**

Fluorescence detection step is at 62°C.

ROX Reference Dye

ROX Reference Dye can be included in the reaction to normalize the fluorescent reporter signal, for instruments that are compatible with that option. ROX is often supplied at a 25µM concentration. Use the following table to determine the amount of ROX to use with a particular instrument:

Instrument	Amount of ROX per 20µl reaction	Final ROX Conc.
ABI 7000, 7300, 7900HT	0.4µl	500 nM
ABI7500 Mx3000P, Mx4000	0.04µl	50 nM

**NOTE**

To accurately pipette 0.04 µl per reaction, dilute ROX 1:10 immediately before use and use 0.4µl of the dilution.

**IMPORTANT**

1. For multiple reactions, prepare a master mix of common components, add the appropriate volume to each tube or plate well on ice, and then add the unique reaction components (*e.g.*, template). Preparation of a master mix is *crucial* in qRT-PCR to reduce pipetting errors.
2. Make sure that all components are at the bottom of the tube/plate; centrifuge briefly if needed.

Data Analyzing

Calibrator and Normalizer

Settling a Calibrator

Before a microRNA relative quantification assay, you should ascertain which sample to be the calibrator. Usually a calibrator is the sample from normal cell or cells without any treatment.

Using U6 snRNA as Normalizer

For relative quantification of microRNA using comparative C_T method, the normalizer is another gene except microRNA, which is often the said housekeeping gene. Subtract the normalizer C_T value from a microRNA C_T value to calculate the ΔC_T for calibrator and samples in each microRNA profiling assay.

IMPORTANT


We recommend normalizer C_T values do not greater than or equal to 30 for any sample. If the normalizer C_T values is consistently equal to or greater than 30, we recommends that you increase the amount of initial total RNA (converted to cDNA) used in each reaction.

Calculating microRNA Expression Ratio Relative to U6 snRNA

Step 1

Designing microRNA relative quantitation assay

It's important to caring two or three repeats for each reaction at least.

Step 2

Run real-time PCR to obtain the C_T value for microRNA and U6 snRNA for each reaction.

Sample 1		Sample 2		Calibrator	
hsa-miR-16	U6 snRNA	hsa-miR-16	U6 snRNA	hsa-miR-16	U6 snRNA
30.40	23.63	24.21	22.66	26.21	24.60
30.35	23.40	24.60	22.56	26.15	24.31
30.41	23.52	24.66	22.48	26.35	24.72

continued on next page

Data Analyzing Continued

IMPORTANT

Attention that the CT value difference between repeat reactions should not greater than 0.5. Otherwise, you should delete the invalid data or carry out the experiment again.

Step 3

Average the C_T values for triplicate wells of the samples and calibrator.

	Sample 1		Sample 2		Calibrator	
	hsa-miR-16	U6 snRNA	hsa-miR-16	U6 snRNA	hsa-miR-16	U6 snRNA
Average C_T	30.39	23.51	24.49	22.56	26.24	24.54

Step 4

Subtract the U6 snRNA values from the miRNAs C_T values to calculate ΔC_T for the calibrator and samples in each reaction.

$$\text{Sample 1 } \Delta C_T = C_T(\text{miR-16}) - C_T(\text{U6snRNA}) = 30.39 - 23.51 = 6.88$$

	Sample 1		Sample 2		Calibrator	
	hsa-miR-16	U6 snRNA	hsa-miR-16	U6 snRNA	hsa-miR-16	U6 snRNA
Average C_T	30.39	23.51	24.49	22.56	26.24	24.54
ΔC_T	6.88	—	1.93	—	1.70	—

Step 5

Subtract the $\Delta C_T(\text{calibrator})$ values from the $\Delta C_T(\text{samples})$ values to calculate their $\Delta\Delta C_T$ (sample) values.

$$\Delta\Delta C_T(\text{sample1}) = \Delta C_T(\text{sample1}) - \Delta C_T(\text{calibrator1}) = 6.88 - 1.70 = 5.18$$

	Sample 1		Sample 2		Calibrator	
	hsa-miR-16	U6 snRNA	hsa-miR-16	U6 snRNA	hsa-miR-16	U6 snRNA
Average C_T	30.39	23.51	24.49	22.56	26.24	24.54
ΔC_T	6.88	—	1.93	—	1.70	—
$\Delta\Delta C_T$	5.18	—	0.23	—	0.00	—

Step 6

Calculate the relative expression ratio.

$$\text{Relative Expression Ratio}_{(\text{sample1})} = 2^{-\Delta\Delta C_T(\text{sample1})} = 2^{-5.18} = 0.027$$

	Sample 1		Sample 2		Calibrator	
	hsa-miR-16	U6 snRNA	hsa-miR-16	U6 snRNA	hsa-miR-16	U6 snRNA
Average C_T	30.39	23.51	24.49	22.56	26.24	24.54
ΔC_T	6.88	—	1.93	—	1.70	—
$\Delta\Delta C_T$	5.18	—	0.23	—	0.00	—
$2^{-\Delta\Delta C_T}$	0.027	—	0.85	—	1.0	—

Appendix

Technical Service

World Wide Web

Visit the GenePharma Web Resource using your World Wide Web browser. At the site, you can:

- ☐ Download manuals in Adobe Acrobat (PDF) format
- ☐ Explore our catalog with full color graphics
- ☐ Get the scoop on our hot new products and special product offers
- ☐ Request catalog and product literature

The Genepharma URL is www.genepharma.com

Quality Control

The product is tested functionally by qRT-PCR using total HeLa RNA as template. Kinetic analysis must demonstrate a linear dose response with decreasing target concentration and detection of U6 snRNA in 1 pg of total HeLa RNA.

Contact Us

For more information or technical assistance, please call, write, fax, or email. Additional international offices are listed on our web page (www.genepharma.com).

Corporate Headquarters:

Genepharma Corporation
1011 Halley Road,
Zhangjiang High-tech Park,
ShangHai, CHINA
Tel: 86-21-51320195
Fax: 86-21-51320295
E-mail: service@genepharma.com

©2005 – 2006 Genepharma Corporation. All rights reserved. For research use only. Not intended for any animal or human therapeutic or diagnostic use.