Lipofectamine® 2000 Reagent

[Learn More](#)
invitrogen™
by *life* technologies™

Package Contents

Catalog Number

- 11668-030 [Buy Now](#)
- 11668-027 [Buy Now](#)
- 11668-019 [Buy Now](#)
- 11668-500 [Buy Now](#)

Size

- 0.3 mL vial
- 0.75 mL vial
- 1.5 mL vial
- 15 mL vial

Storage Conditions

Store at 4°C (do not freeze).

Required Materials

- Plasmid DNA (0.5–5 µg/µL stock)
- Opti-MEM® Reduced Serum Medium
- Eppendorf tubes

Timing

Preparation: 10 minutes
Incubation: 5 minutes
Final Incubation: 1–3 days

Selection Guide

[Lipofectamine® Reagents](#)

Go online to view related products.

Product Description

- Lipofectamine® 2000 Reagent is a proprietary formulation for transfecting nucleic acids into a wide range of eukaryotic cells.

Important Guidelines

- DNA-Lipofectamine® 2000 complexes must be made in serum-free medium such as Opti-MEM® Reduced Serum Medium and can be added directly to cells in culture medium, in the presence or absence of serum/antibiotic.
- It is not necessary to remove complexes or change/add medium after transfection.
- The amount of Lipofectamine® 2000 Reagent required for successful transfection varies depending on the cell type and passage number. Start any new transfection by testing the recommended four concentrations of Lipofectamine® 2000 Reagent to determine an optimum amount.

Online Resources

Visit our [product page](#) for additional information and protocols. For support, visit www.lifetechnologies.com/support.

Protocol Outline

- Plate cells so they will be 70–90% confluent at the time of transfection.
- Prepare plasmid DNA-lipid complexes.
- Add DNA-lipid complexes to cells.

Lipofectamine® 2000 DNA Transfection Reagent Protocol

- i** See page 2 to view a typical DNA transfection procedure.

Component	96-well	24-well	6-well
Final DNA per well	100 ng	500 ng	2500 ng
Final Lipofectamine® 2000 Reagent per well	0.2–0.5 µL	1.0–2.5 µL	5.0–12.5 µL

Co-Transfection of Plasmid DNA and siRNA

Transfect plasmid DNA and siRNA at the same time using Lipofectamine® 2000 Reagent by adding 30 pmol (~0.6 µg) of siRNA per 1 µg of DNA.

mRNA Transfection

mRNA can be transfected in a 24-well plate using Lipofectamine® 2000 Reagent by adding 0.5–1 µg of mRNA per well.

i Photograph of Expected Results

i Scaling Up or Down Transfections

i Limited Product Warranty and Disclaimer Details

Lipofectamine® 2000 DNA Transfection Reagent Protocol

Transfect cells according to the following chart. Volumes are given on a per-well basis. **Each reaction mix is sufficient for triplicate (96-well), duplicate (24-well), and single well (6-well) transfections, and accounts for pipetting variations.** Adjust the amounts of components according to your tissue culture format. For additional information on scaling your transfection reaction, see page 1.

Timeline		Steps	Procedure Details			
Day 0	1	Seed cells to be 70–90% confluent at transfection	Component	96-well	24-well	6-well
	2	Dilute four amounts of Lipofectamine® Reagent in Opti-MEM® Medium	Adherent cells	$1-4 \times 10^4$	$0.5-2 \times 10^5$	$0.25-1 \times 10^6$
	3	Dilute DNA in Opti-MEM® Medium	Opti-MEM® Medium	25 μ L $\times$ 4	50 μ L $\times$ 4	150 μ L $\times$ 4
Day 1	4	Add diluted DNA to diluted Lipofectamine® 2000 Reagent (1:1 ratio)	Lipofectamine® 2000 Reagent	i 1, 1.5, 2, 2.5 μ L	i 2, 3, 4, 5 μ L	i 6, 9, 12, 15 μ L
	5	Incubate	Opti-MEM® Medium	125 μ L	250 μ L	700 μ L
	6	Add DNA-lipid complex to cells	DNA (0.5–5 μ g/ μ L)	2.5 μ g	5 μ g	14 μ g
	7	Visualize/analyze transfected cells	Diluted DNA Total	25 μ L	50 μ L	150 μ L
Day 2–4			Diluted Lipofectamine® 2000 Reagent	25 μ L	50 μ L	150 μ L
			Incubate for 5 minutes at room temperature.			
			Component	96-well	24-well	6-well
			DNA-lipid complex per well	10 μ L	50 μ L	250 μ L
			Final DNA used per well	100 ng	500 ng	2500 ng
			Final Lipofectamine® 2000 Reagent used per well	0.2–0.5 μ L	1.0–2.5 μ L	5.0–12.5 μ L
			Incubate cells for 1–3 days at 37°C. Then analyze transfected cells.			

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